

GPR101 Antibody, Biotin conjugated



PACO39341

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Orphan receptor.

Gene ID:

GPR101

Uniprot

Q96P66

Synonyms:

Probable G-protein coupled receptor 101, GPR101

Immunogen:

Recombinant Human Probable G-protein coupled receptor 101 protein (218-399AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A