

PGLYRP2 Antibody, Biotin conjugated



PACO39345

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a scavenger role by digesting biologically active peptidoglycan (PGN) into biologically inactive fragments. Has no direct bacteriolytic activity.

Gene ID:

PGLYRP2

Uniprot

Q96PD5

Synonyms:

N-acetylmuramoyl-L-alanine amidase (EC 3.5.1.28) (Peptidoglycan recognition protein 2) (Peptidoglycan recognition protein long) (PGRP-L), PGLYRP2, PGLYRPL PGRPL

Immunogen:

Recombinant Human N-acetylmuramoyl-L-alanine amidase protein (22-576AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A