

CA13 Antibody, FITC conjugated



PACO39388

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Reversible hydration of carbon dioxide.

Gene ID:

CA13

Uniprot

Q8N1Q1

Synonyms:

Carbonic anhydrase 13 (EC 4.2.1.1) (Carbonate dehydratase XIII) (Carbonic anhydrase XIII) (CA-XIII), CA13

Immunogen:

Recombinant Human Carbonic anhydrase 13 protein (1-262AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A