

TMEM256 Antibody



PACO39394

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

extracellular exosome.

Gene ID:

TMEM256

Uniprot

Q8N2U0

Synonyms:

Transmembrane protein 256, TMEM256, C17orf61


Immunogen:

Recombinant Human Transmembrane protein 256 protein (30-63AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of MCF-7 cells using PACO39394 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).