

TCP11L2 Antibody, Biotin conjugated



PACO39417

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TCP11L2

Uniprot

Q8N4U5

Synonyms:

T-complex protein 11-like protein 2, TCP11L2

Immunogen:

Recombinant Human T-complex protein 11-like protein 2 protein (1-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A