

SAMD12 Antibody



PACO39434

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

SAMD12

Uniprot

Q8N8I0

Synonyms:

Sterile alpha motif domain-containing protein 12 (SAM domain-containing protein 12), SAMD12


Immunogen:

Recombinant Human Sterile α motif domain-containing protein 12 protein (1-201AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39434 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO39434 at dilution of 1:100.