

ICA1L Antibody



PACO39446

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

ICA1L

Uniprot

Q8NDH6

Synonyms:

Islet cell autoantigen 1-like protein (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 14 protein) (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 15 protein), ICA1L, ALS2CR14 ALS2CR15


Immunogen:

Recombinant Human Islet cell autoantigen 1-like protein (1-187AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: ICA1L antibody at 1.2 μ g/ml + MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 55, 22 kDa. Observed band size: 55 kDa.