

NUP37 Antibody



PACO39454

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the Nup107-160 subcomplex of the nuclear pore complex (NPC). The Nup107-160 subcomplex is required for the assembly of a functional NPC. The Nup107-160 subcomplex is also required for normal kinetochore microtubule attachment, mitotic progression and chromosome segregation.

Gene ID:

NUP37

Uniprot

Q8NFH4

Synonyms:

Nucleoporin Nup37 (p37) (Nup107-160 subcomplex subunit Nup37), NUP37

Immunogen:

Recombinant Human Nucleoporin Nup37 protein (1-326AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A