

# VN1R2 Antibody, HRP conjugated



PACO39459

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Putative pheromone receptor.

**Gene ID:**

VN1R2

**Uniprot**

Q8NFZ6

**Synonyms:**

Vomer nasal type-1 receptor 2 (G-protein coupled receptor GPCR25) (hGPCR25) (V1r-like receptor 2), VN1R2, V1RL2

**Immunogen:**

Recombinant Human Vomer nasal type-1 receptor 2 protein (115-170AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A