

VN1R2 Antibody, Biotin conjugated



PACO39461

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Putative pheromone receptor.

Gene ID:

VN1R2

Uniprot

Q8NFZ6

Synonyms:

Vomer nasal type-1 receptor 2 (G-protein coupled receptor GPCR25) (hGPCR25) (V1r-like receptor 2), VN1R2, V1RL2

Immunogen:

Recombinant Human Vomer nasal type-1 receptor 2 protein (115-170AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A