

# CNIH3 Antibody, Biotin conjugated



PACO39469

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Regulates the trafficking and gating properties of AMPA-selective glutamate receptors (AMPA receptors). Promotes their targeting to the cell membrane and synapses and modulates their gating properties by regulating their rates of activation, deactivation and desensitization.

**Gene ID:**

CNIH3

**Uniprot**

Q8TBE1

**Synonyms:**

Protein cornichon homolog 3 (CNIH-3) (Cornichon family AMPA receptor auxiliary protein 3), CNIH3

**Immunogen:**

Recombinant Human Protein cornichon homolog 3 protein (94-138AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A