

PIGX Antibody, FITC conjugated



PACO39472

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Essential component of glycosylphosphatidylinositol-mannosyltransferase 1 which transfers the first of the 4 mannoses in the GPI-anchor precursors during GPI-anchor biosynthesis. Probably acts by stabilizing the mannosyltransferase PIGM.

Gene ID:

PIGX

Uniprot

Q8TBF5

Synonyms:

Phosphatidylinositol-glycan biosynthesis class X protein (PIG-X), PIGX

Immunogen:

Recombinant Human Phosphatidylinositol-glycan biosynthesis class X protein (42-230AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A