

C16orf78 Antibody




PACO39494

Product Information

Size:	Protein Background:
50ug	nucleus.
Reactivity:	Gene ID:
Human	C16orf78
Source:	Uniprot
Rabbit	Q8WTQ4
Isotype:	Synonyms:
IgG	Uncharacterized protein C16orf78, C16orf78
Applications:	Immunogen:
ELISA, IHC	Recombinant Human Uncharacterized protein C16orf78 protein (1-265AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39494 at dilution of 1:100.