

TBC1D22A Antibody



PACO39498

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

May act as a GTPase-activating protein for Rab family protein(s).

Gene ID:

TBC1D22A

Uniprot

Q8WUA7

Synonyms:

TBC1 domain family member 22A, TBC1D22A, C22orf4


Immunogen:

Recombinant Human TBC1 domain family member 22A protein (1-280AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TBC1D22A antibody at 18 μ g/ml. Lane 1: Colo320 whole cell lysate. Lane 2: HepG2 whole cell lysate. Lane 3: Rat skeletal muscle tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 60, 51, 37, 54 kDa. Observed band size: 60 kDa.