

TBC1D22A Antibody, HRP conjugated



PACO39499

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May act as a GTPase-activating protein for Rab family protein(s).

Gene ID:

TBC1D22A

Uniprot

Q8WUA7

Synonyms:

TBC1 domain family member 22A, TBC1D22A, C22orf4

Immunogen:

Recombinant Human TBC1 domain family member 22A protein (1-280AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A