

GTSF1 Antibody, FITC conjugated



PACO39508

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GTSF1

Uniprot

Q8WW33

Synonyms:

Gametocyte-specific factor 1 (Protein FAM112B), GTSF1, FAM112B

Immunogen:

Recombinant Human Gametocyte-specific factor 1 protein (1-167AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A