

RTN4IP1 Antibody



PACO39510

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Appears to be a potent inhibitor of regeneration following spinal cord injury.

Gene ID:

RTN4IP1

Uniprot

Q8WWV3

Synonyms:

Reticulon-4-interacting protein 1, mitochondrial (NOGO-interacting mitochondrial protein), RTN4IP1, NIMP


Immunogen:

Recombinant Human Reticulon-4-interacting protein 1, mitochondrial protein (41-396AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human appendix tissue using PACO39510 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39510 at dilution of 1:100.