

# GALM Antibody, FITC conjugated



PACO39556

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Mutarotase converts alpha-aldose to the beta-anomer. It is active on D-glucose, L-arabinose, D-xylose, D-galactose, maltose and lactose.

**Gene ID:**

GALM

**Uniprot**

Q96C23

**Synonyms:**

Aldose 1-epimerase (EC 5.1.3.3) (Galactose mutarotase), GALM

**Immunogen:**

Recombinant Human Aldose 1-epimerase protein (2-342AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A