

SHD Antibody



PACO39598

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May function as an adapter protein.

Gene ID:

SHD

Uniprot

Q96IW2

Synonyms:

SH2 domain-containing adapter protein D, SHD


Immunogen:

Recombinant Human SH2 domain-containing adapter protein D protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO39598 at dilution of 1:100.