

RBM41 Antibody



PACO39602

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

May bind RNA.

Gene ID:

RBM41

Uniprot

Q96IZ5

Synonyms:

RNA-binding protein 41 (RNA-binding motif protein 41), RBM41


Immunogen:

Recombinant Human RNA-binding protein 41 protein (1-413AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: RBM41 antibody at 1.5 μ g/ml + HeLa whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 48, 47, 36 kDa. Observed band size: 48 kDa.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO39602 at dilution of 1:100.