

MRGPRX3 Antibody, Biotin conjugated



PACO39613

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Orphan receptor. Probably involved in the function of nociceptive neurons. May regulate nociceptor function and/or development, including the sensation or modulation of pain. Potently activated by enkephalins.

Gene ID:

MRGPRX3

Uniprot

Q96LB0

Synonyms:

Mas-related G-protein coupled receptor member X3 (Sensory neuron-specific G-protein coupled receptor 1/2), MRGPRX3, MRGX3 SNSR1 SNSR2

Immunogen:

Recombinant Human Mas-related G-protein coupled receptor member X3 protein (276-322AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A