

AMDHD1 Antibody, FITC conjugated



PACO39616

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

AMDHD1

Uniprot

Q96NU7

Synonyms:

Probable imidazolonepropionase (EC 3.5.2.7) (Amidohydrolase domain-containing protein 1), AMDHD1

Immunogen:

Recombinant Human Probable imidazolonepropionase protein (1-426AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A