

KLHL34 Antibody



PACO39634

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Cul3-RING ubiquitin ligase complex, extracellular space, protein ubiquitination.

Gene ID:

KLHL34

Uniprot

Q8N239

Synonyms:

Kelch-like protein 34, KLHL34


Immunogen:

Recombinant Human Kelch-like protein 34 protein (401-644AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO39634 at dilution of 1:100.