

PACRGL Antibody



PACO39658

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PACRGL

Uniprot

Q8N7B6

Synonyms:

PACRG-like protein, PACRGL, C4orf28

Immunogen:

Recombinant Human PACRG-like protein (1-130AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A