

PRELID2 Antibody



PACO39670

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

mitochondrial intermembrane space, phosphatidic acid, transporter activity, phospholipid transport.

Gene ID:

PRELID2

Uniprot

Q8N945

Synonyms:

PRELI domain-containing protein 2, PRELID2

Immunogen:

Recombinant Human PRELI domain-containing protein 2 protein (1-148AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A