

PACO39674

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role into cellular response to oxidative stress. Decreases cell proliferation.

Gene ID:

EIF1AD

Uniprot

Q8N9N8

Synonyms:

Probable RNA-binding protein EIF1AD (Eukaryotic translation initiation factor 1A domain-containing protein) (Haponin), EIF1AD

Immunogen:

Recombinant Human Probable RNA-binding protein EIF1AD protein (1-165AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A