

NHLRC2 Antibody



PACO39686

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

Gene ID:

NHLRC2

Uniprot

Q8NBF2

Synonyms:

NHL repeat-containing protein 2, NHLRC2


Immunogen:

Recombinant Human NHL repeat-containing protein 2 protein (1-367AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: NHLRC2 antibody at 8 μ g/ml + U251 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 80, 40 kDa. Observed band size: 80 kDa.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO39686 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO39686 at dilution of 1:100.