

# TTC16 Antibody



PACO39694

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

TTC16

**Uniprot**

Q8NEE8

**Synonyms:**

Tetratricopeptide repeat protein 16 (TPR repeat protein 16), TTC16

**Immunogen:**


Recombinant Human Tetratricopeptide repeat protein 16 protein (551-873AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39694 at dilution of 1:100.