

CFAP61 Antibody



PACO39710

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

CFAP61

Uniprot

Q8NHU2

Synonyms:

Cilia- and flagella-associated protein 61, CFAP61, C20orf26


Immunogen:

Recombinant Human Cilia- and flagella-associated protein 61 protein (1-200AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39710 at dilution of 1:100.



Immunohistochemistry analysis of human pancreatic cancer using PACO39710 at dilution of 1:100.



Immunohistochemistry analysis of human liver cancer using PACO39710 at dilution of 1:100.