

# FCGRT Antibody, Biotin conjugated



PACO39717

---

## Product Information

**Size:**

50ug

**Reactivity:**

Macaca fascicularis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds to the Fc region of monomeric immunoglobulins gamma. Mediates the uptake of IgG from milk. Possible role in transfer of immunoglobulin G from mother to fetus.

**Gene ID:**

FCGRT

**Uniprot**

Q8SPV9

**Synonyms:**

IgG receptor FcRn large subunit p51 (FcRn) (IgG Fc fragment receptor transporter alpha chain) (Neonatal Fc receptor), FCGRT, FCRN

**Immunogen:**

Recombinant Macaca fascicularis IgG receptor FcRn large subunit p51 protein (24-297AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A