

FIBIN Antibody



PACO39718

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FIBIN

Uniprot

Q8TAL6

Synonyms:

Fin bud initiation factor homolog, FIBIN

Immunogen:

Recombinant Human Fin bud initiation factor homolog protein (19-211AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A