

PACO39722

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Non-mitochondrial 3-hydroxymethyl-3-methylglutaryl-CoA lyase that catalyzes a cation-dependent cleavage of (S)-3-hydroxy-3-methylglutaryl-CoA into acetyl-CoA and acetoacetate, a key step in ketogenesis, the products of which support energy production in nonhepatic animal tissues.

Gene ID:

HMGCLL1

Uniprot

Q8TB92

Synonyms:

3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic (EC 4.1.3.4) (3-hydroxymethyl-3-methylglutaryl-CoA lyase-like protein 1) (Endoplasmic reticulum 3-hydroxymethyl-3-methylglutaryl-CoA lyase) (er-cHL), HMGCLL1

Immunogen:

Recombinant Human 3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic protein (2-340AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A