

FAM174A Antibody



PACO39726

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

FAM174A

Uniprot

Q8TBP5

Synonyms:

Membrane protein FAM174A (Hepatitis C virus NS5A-transactivated protein 6) (HCV NS5A-transactivated protein 6) (Transmembrane protein 157), FAM174A, NS5ATP6, TMEM157


Immunogen:

Recombinant Human Membrane protein FAM174A protein (41-123AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO39726 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39726 at dilution of 1:100.