

TEX47 Antibody



PACO39730

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

TEX47

Uniprot

Q8TBZ9

Synonyms:

Testis-expressed protein 47, TEX47, C7orf62


Immunogen:

Recombinant Human Testis-expressed protein 47 protein (1-253AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TEX47 antibody at 8 μ g/ml. Lane 1: Jurkat whole cell lysate. Lane 2: A549 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 30 kDa. Observed band size: 30 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39730 at dilution of 1:100.