

Plet1 Antibody, Biotin conjugated



PACO39741

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Modulates leading keratinocyte migration and cellular adhesion to matrix proteins during a wound-healing response and promotes wound repair. May play a role during trichilemmal differentiation of the hair follicle.

Gene ID:

Plet1

Uniprot

Q8VEN2

Synonyms:

Placenta-expressed transcript 1 protein (Antigen mAgK114), Plet1

Immunogen:

Recombinant Mouse Placenta-expressed transcript 1 protein (28-218AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A