

FAM222B Antibody



PACO39742


Product Information

Size:	Protein Background:
50ug	nucleoplasm.
Reactivity:	Gene ID:
Human, Mouse	FAM222B
Source:	Uniprot
Rabbit	Q8WU58
Isotype:	Synonyms:
IgG	Protein FAM222B, FAM222B, C17orf63
Applications:	Immunogen:
ELISA, WB, IHC	Recombinant Human Protein FAM222B protein (351-562AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:1000-1:5000, IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: FAM222B antibody at 2 μ g/ml + Mouse liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 60 kDa. Observed band size: 60 kDa.



Immunohistochemistry of paraffin-embedded human bladder cancer using PACO39742 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon tissue using PACO39742 at dilution of 1:100.