

CPLX3 Antibody, Biotin conjugated



PACO39757

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Positively regulates a late step in synaptic vesicle exocytosis.

Gene ID:

CPLX3

Uniprot

Q8WVH0

Synonyms:

Complexin-3 (Complexin III) (CPX III), CPLX3

Immunogen:

Recombinant Human Complexin-3 protein (1-155AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A