

sodB Antibody



PACO39770

Product Information

Size:

50ug

Reactivity:

Neosartorya fumigata

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Destroys superoxide anion radicals which are normally produced within the cells and which are toxic to biological systems.

Gene ID:

sodB

Uniprot

Q92450

Synonyms:

Superoxide dismutase [Mn], mitochondrial (EC 1.15.1.1) (allergen Asp f 6), sodB

Immunogen:

Recombinant Neosartorya fumigata Superoxide dismutase [Mn], mitochondrial protein (1-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A