

GTF3C6 Antibody, FITC conjugated



PACO39784

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in RNA polymerase III-mediated transcription. Integral, tightly associated component of the DNA-binding TFIIIC2 subcomplex that directly binds tRNA and virus-associated RNA promoters.

Gene ID:

GTF3C6

Uniprot

Q969F1

Synonyms:

General transcription factor 3C polypeptide 6 (Transcription factor IIIC 35 kDa subunit) (TFIIIC 35 kDa subunit) (TFIIIC35) (Transcription factor IIIC subunit 6), GTF3C6, C6orf51

Immunogen:

Recombinant Human General transcription factor 3C polypeptide 6 protein (2-213AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A