

# TMBIM1 Antibody, FITC conjugated



PACO39792

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Negatively regulates aortic matrix metalloproteinase-9 (MMP9) production and may play a protective role in vascular remodeling.

**Gene ID:**

TMBIM1

**Uniprot**

Q969X1

**Synonyms:**

Protein lifeguard 3 (Protein RECS1 homolog) (Transmembrane BAX inhibitor motif-containing protein 1), TMBIM1, LFG3 RECS1

**Immunogen:**

Recombinant Human Protein lifeguard 3 protein (1-109AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A