

MRPL53 Antibody, HRP conjugated



PACO39819

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MRPL53

Uniprot

Q96EL3

Synonyms:

39S ribosomal protein L53, mitochondrial (L53mt) (MRP-L53) (Mitochondrial large ribosomal subunit protein mL53), MRPL53

Immunogen:

Recombinant Human 39S ribosomal protein L53, mitochondrial protein (8-113AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A