

S100A16 Antibody, Biotin conjugated



PACO39829

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Calcium-binding protein. Binds one calcium ion per monomer.

Gene ID:

S100A16

Uniprot

Q96FQ6

Synonyms:

Protein S100-A16 (Aging-associated gene 13 protein) (Protein S100-F) (S100 calcium-binding protein A16), S100A16, S100F

Immunogen:

Recombinant Human Protein S100-A16 protein (1-103AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A