

FAM163A Antibody



PACO39842

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

FAM163A

Uniprot

Q96GL9

Synonyms:

Protein FAM163A (Cebelin) (Neuroblastoma-derived secretory protein), FAM163A, C1orf76 NDSP


Immunogen:

Recombinant Human Protein FAM163A protein (27-167AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39842 at dilution of 1:100.