

MAGEB10 Antibody



PACO39866

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:**Gene ID:**

MAGEB10

Uniprot

Q96LZ2

Synonyms:

Melanoma-associated antigen B10 (MAGE-B10 antigen), MAGEB10

Immunogen:

Recombinant Human Melanoma-associated antigen B10 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: MAGEB10 antibody at 12 μ g/ml + Mouse heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 39 kDa. Observed band size: 39, 35 kDa.