

PRRC1 Antibody



PACO39870

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

PRRC1

Uniprot

Q96M27

Synonyms:

Protein PRRC1 (Proline-rich and coiled-coil-containing protein 1), PRRC1


Immunogen:

Recombinant Human Protein PRRC1 protein (200-445AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO39870 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO39870 at dilution of 1:100.