

SFXN2 Antibody, FITC conjugated



PACO39880

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Potential iron transporter.

Gene ID:

SFXN2

Uniprot

Q96NB2

Synonyms:

Sideroflexin-2, SFXN2

Immunogen:

Recombinant Human Sideroflexin-2 protein (1-99AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A