

TNFRSF13C Antibody, Biotin conjugated



PACO39889

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

B-cell receptor specific for TNFSF13B/TALL1/BAFF/BLyS. Promotes the survival of mature B-cells and the B-cell response.

Gene ID:

TNFRSF13C

Uniprot

Q96RJ3

Synonyms:

Tumor necrosis factor receptor superfamily member 13C (B-cell-activating factor receptor) (BAFF receptor) (BAFF-R) (BLyS receptor 3) (CD antigen CD268), TNFRSF13C, BAFFR BR3

Immunogen:

Recombinant Human Tumor necrosis factor receptor superfamily member 13C protein (1-78AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A