

# ARHGEF39 Antibody



PACO39918

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Promotes cell proliferation.

**Gene ID:**

ARHGEF39

**Uniprot**

Q8N4T4

**Synonyms:**

Rho guanine nucleotide exchange factor 39, ARHGEF39, C9orf100

**Immunogen:**

Recombinant Human Rho guanine nucleotide exchange factor 39 protein (22-197AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A