

FAM89B Antibody, FITC conjugated



PACO39928

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Negatively regulates TGF-beta-induced signaling; in cooperation with SKI prevents the translocation of SMAD2 from the nucleus to the cytoplasm in response to TGF-beta.

Gene ID:

FAM89B

Uniprot

Q8N5H3

Synonyms:

Leucine repeat adapter protein 25, FAM89B, Lrap25

Immunogen:

Recombinant Human Leucine repeat adapter protein 25 protein (1-162AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A