

ACBD7 Antibody



PACO39934

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds medium- and long-chain acyl-CoA esters.

Gene ID:

ACBD7

Uniprot

Q8N6N7

Synonyms:

Acyl-CoA-binding domain-containing protein 7, ACBD7

Immunogen:

Recombinant Human Acyl-CoA-binding domain-containing protein 7 protein (1-88AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A